

Incert[®]-S

For the prophylaxis of peridural adhesions

"Ease approach in case of revision surgery."

Peridural scarring is a natural consequence of most spine surgery. Though difficult to prove in humans, a number of studies have suggested a link between post-operative peridural scarring and radicular pain.

The presence of peridural adhesions may also compromise success of fusion procedures, and cause extraordinary surgical challenges in the event that secondary operations are required at the same site later in life.

For this reason, surgeons have sought a consistent, effective method which complements current and emerging surgical techniques to prevent post-operative scarring.

Incert[®]-S is a clear anti-adhesion viscous barrier gel, based on a patented biomaterial, consisting of highly purified, chemically modified, high molecular weight Hyaluronic Acid (HA).

Incert[®]-S provides the surgeon with a biocompatible, bioabsorbable anti-adhesion barrier gel to be used as an adjunct to surgery during laminectomy, laminotomy and discectomy procedures.

Incert[®]-S coats the sensitive neural structures during a short period after surgical closure and provides protection from the formation of peridural fibrosis and dural adhesions that may otherwise contribute to post-operative radicular pain and/or neurological dysfunction during the formation period.

Incert[®]-S comes pre-filled in a syringe, with an applicator, to facilitate easy application to the spinal neural structures usually accessed during laminectomy and discectomy procedures.

Incert[®]-S offers the surgeon confidence through its:

- Biocompatibility** ensuring the patients safety
- Bioabsorbable** leaving no residual foreign elements to impact surgical intervention
- Adherence** allowing it to remain in place for the desired effect
- Ease of use**

Ordering Information

Product Code: Incert-S Description: 2ml Syringe with applicator

*Data on file

Authorized Distributor (Malaysia) :



CIMed Healthcare Sdn Bhd

(889001-X)

Enrich life through Innovation

2-112, Jalan Prima SG 3/2, Taman Prima Sri Gombak,
68100 Batu Caves, Selangor Darul Ehsan, Malaysia.
Tel : +603-61888 089 Fax : +603-61888 090
Email : info@cimedhealthcare.com www.cimedhealthcare.com



Manufactured by Anika Therapeutics, Woburn, MA.

